



## Full Automatic Gas Pycnometer: Macpycno<sup>®</sup>

Full automatic gas pycnometer for reliable, precise and quick analysis of true volume and real density of solids by using the gas displacement method.



The Macpycno<sup>®</sup> provides reliable and precise results in a short amount of time. Additionally, with its high repeatability it is most suitable for laboratories and quality control.

To set up the measurement of a sample there are only a few steps needed. The actual measurement itself is fully automatic. Macpycno<sup>®</sup> is also equipped with Peltier device to monitor and adjust the temperature of the chambers.


A high range of various samples can be analyzed from nano materials to components and many others.

### Measurement principle

The measurement principle relies on the well known gas displacement method to measure the true volume and density of solids. The Macpycno<sup>®</sup> uses helium as the displacement gas due to its size and generally inert behavior.

#### Principle of Operation

- (1) Sample chamber (empty) and reference chamber are purged from atmospheric air with helium.
- (2) Valve 2, 3 are closed and the sample is put inside the sample chamber. The pressure in the sample chamber will rise to a certain level.
- (3) If the max. pressure is reached, valve 1 is closed and valve 2 is opened. A pressure equalization between sample and reference chamber will happen, that is used to calculate the true volume of the sample.



### Specifications

Principle	Gas displacement method
Manufacturer's country	Japan
Displacement gas	Helium
Repeatability	± 0.03%
Measurement Cell Volume	Sample pod : 20cc, 40cc, 60cc ※Option : 3cc, 110cc
Calibration	Stainless steel ball for calibration
Interface	Color touch panel
Chamber temperature control	18~30°C (automatically controlled through Peltier device) Measurement pressure self control (Option)
Size	W450 × D410 × H250
Weight	26kg
Note	Possible to output the measurement data to USB stick as CSV file Possible to manage and analyze the data by PC.